



CERTIFICATE OF COMPARISON

NAME

CERT NO.

GC393

PRODUCT DESCRIPTION

S12 STAGE MICROMETER SCALE

DATE ISSUED

11 March 2015

SERIAL NUMBER

PM-ID00014

MEASURED INTERVALS

0.02mm intervals

The above product has been measured in our laboratories and compared with a NAMAS certified scale as described below. The accumulative error over the total length is within **±0.0005mm.**

EXAMINED BY *STAN KOWAL*

APPROVED

UNCERTAINTY OF MEASUREMENT

+/- 1 micron

Calibration Equipment:

Tetrascan

Calibration Standard used:

S12 Stage Micrometer Scale

Serial No:

A387

Certificate No:

08A032/95019/ML72/41 (Calibrated by NPL)

Telephone: +44 (0)1732 864111

Facsimile: +44 (0)1732 865544