

Specimen Stubs for mounting calibration standards

Calibration standards can be either purchased unmounted or mounted on a specific specimen stub. The mounting stubs offered are:

-1	#10-002012	Standard Pin Stub Ø12.7mm, 9.5mm pinH
-2	#10-003012	Standard Pin Stub Ø12.7mm, short 6mm pinH
-3	#10-002025	Large Pin Stub Ø25mm, 9.5mm pinH
-4	#10-003025	Large Pin Stub Ø25mm, 6mm pinH
-5	#10-005110	JEOL Ø9.5x9.5mm Cylinder Stub
-6	#10-005112	JEOL Ø12.2x10mm Cylinder Stub
-7	#10-005125	JEOL Ø25x10mm Cylinder Stub
-8	#10-004015	Hitachi Ø15x6mmxM4 Cylinder Stub
-9	#10-004025	Hitachi Ø25x6mmxM4 Cylinder Stub
-10	you supply	Custom Stub mounting, stub supplied by customer
-11	#51-001001	Black metal slide

If you wish to use a calibration standard across different SEM platforms, consider using the EM-Tec SEM stub adapters.

If you wish to have the calibration standard mounted on a different SEM stub, you can choose the custom stub mounting. You have to supply the mount for custom mounting.

TSB 062115-2 Specimen stubs for mounting calibration standards

TSB 062115-1 Specimen stubs for mounting calibration standards 2015-06-21